

## Evita 4

### Technical Service Bulletin # 36

|                          |   |               |  |
|--------------------------|---|---------------|--|
| <b>Re:</b>               | <b>New software version 3.21</b>  |               |  |
| <b>Update:</b>           | May 22, 2000  |               |  |
| <b>Reference Doc:</b>    | IDM nos. 36 / 09.99, 35 / 08.99, 32 / 05.99, and 28 / 02.99   |               |  |
| <b>Reason:</b>           | <p>Product improvement. This software version 3.21 replaces all older versions in the USA and Canada.</p> <p>These are the most important changes for the operator if the device is upgraded from software version 2.n to 3.21:</p> <ul style="list-style-type: none"><li>• New "Night &amp; Day" key to dim the display screen.</li><li>• Additional mode ILV (independent lung ventilation).</li><li>• New "Exp. Hold" key to prolong the expiration phase.</li><li>• The nebulization function is now available in the pediatric mode.</li><li>• Oxygenation for bronchial suction in the pediatric mode, increase of O<sub>2</sub> concentration by 25 % relative (previously at 100 Vol. %).</li></ul> |               |  |
| <b>Solution:</b>         | Installation of new software version SW 3.21.   |               |  |
| <b>Devices affected:</b> | <p>Recommended for all Evita 4 units with software version 2.n or 1.n.</p> <p>Mandatory on units with installed software version 3.05, 3.06, and 3.10.</p> <p>These software versions were not delivered to the USA and Canada.</p>   |               |  |
| <b>When:</b>             | When requested by Customer, next Service Call or Preventive Maintenance.  |               |  |
| <b>Additional Info:</b>  | <p>After downloading the new software, perform a device check in stand-by. The message "Loss of data" may occur shortly after first power on after this download and is not a device malfunction. The message will disappear after the internal data backup is completed.</p> <p>The new functions of software version 3.21 will require a clinical inservice if the previous version was software version 1.n or 2.n.</p>  |               |  |
| <b>Cost:</b>             | Covered by Customer.  |               |  |
| <b>Ordering Info:</b>    | Evita 4 software version 3.21 upgrade kit   | P/N 84 14 589 |  |

The kit consists of Operating Instructions, option label, disk with SW 3.21, and new keypad labels.

**Distribution:** Dräger Service Personnel and Authorized Service Organizations for CCS products.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-721-5789

Dräger Medical, Inc.  
Technical Product Manager